

System Name	Size	Product Name	Use with main product as combination (Product or Accessory)					Instrument Set code	Instrument Set Name	Instrument Set Details	Mention if used on joints (Implants used on joints will affect the adjacent body part)	Name of surgical technique used
			Screws (If Applicable) With Quantity	Pins (If Applicable) With Quantity	Wire (If Applicable) With Quantity	Bolt /END CAP (If Applicable) With Quantity	Fixators (If Applicable) With Quantity					
TIBIA NAIL LOCKING SYSTEM	8 mm x (280 mm – 380 mm)With Diff. of 20 mm	INTERLOCKING NAIL FOR TIBIA	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt3.9 mm (4 Pcs).& Standard End cap for 8 mm & 9 mm	-	287.31	Tibia/Femur Set	As per Annexure A	-	Interlocking nail system YIPL/NL/05
	9 mm x (280 mm – 380 mm)With Diff. of 20 mm	INTERLOCKING NAIL FOR TIBIA	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4pcs)& Standard End cap for 8 mm & 9 mm	-	287.31	Tibia/Femur Set	As per Annexure A	-	Interlocking nail system YIPL/NL/05
	10 mm x (280 mm – 380 mm)With Diff. of 20 mm	INTERLOCKING NAIL FOR TIBIA	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4pcs)& Standard End cap for 10 mm & 11 mm	-	287.31	Tibia/Femur Set	As per Annexure A	-	Interlocking nail system YIPL/NL/05

DISTAL TIBIA NAIL LOCKING SYSTEM	8 mm x (280 mm – 380 mm)With Diff. of 20 mm	Distal Tibia Nail - Cannulated	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt3.9 mm (4 Pcs).	-	287.31	Tibia/Femur Set	As per Annexure B	-	Interlocking nail system YIPL/NL/05
	9 mm x (280 mm – 380 mm)& Diff. of 20 mm	Distal Tibia Nail - Cannulated	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 pcs)	-	287.31	Tibia/Femur Set	As per Annexure B	-	Interlocking nail system YIPL/NL/05
	10 mm x (280 mm – 380 mm)With Diff. of 20 mm	Distal Tibia Nail - Cannulated	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 pcs)	-	287.31	Tibia/Femur Set	As per Annexure B	-	Interlocking nail system YIPL/NL/05
EXCEL NAIL LOCKING SYSTEM	8 mm x (280 mm – 380 mm)With Diff. of 20 mm	EXCEL NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt3.9 mm (4 Pcs).	-	289.25	G.T.N./EXCEL SET	As per Annexure C	-	Excel nail YIPL/NL/03
	9 mm x (280 mm – 380 mm)With Diff. of 20 mm	EXCEL NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4pcs)	-	289.25	G.T.N./EXCEL SET	As per Annexure C	-	Excel nail YIPL/NL/03

	10 mm x (280 mm – 380 mm)With Diff. of 20 mm	EXCEL NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4pcs)	-	289.25	G.T.N./ EXCEL SET	As per Annexure C	-	Excel nail YIPL/NL/ 03
EXCEL - 1NAIL LOCKING SYSTEM	8 mm x (280 mm – 380 mm)With Diff. of 20 mm	EXCEL 1 NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs)	-	289.25	G.T.N./ EXCEL SET	As per Annexure C	-	Excel nail YIPL/NL/ 03
	9 mm x (280 mm – 380 mm)With Diff. of 20 mm	EXCEL 1 NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs)	-	289.25	G.T.N./ EXCEL SET	As per Annexure C	-	Excel nail YIPL/NL/ 03
	10 mm x (280 mm – 380 mm)With Diff. of 20 mm	EXCEL 1 NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs)	-	289.25	G.T.N./ EXCEL SET	As per Annexure C	-	Excel nail YIPL/NL/ 03
EXCEL D NAIL LOCKING SYSTEM	8 mm x (280 mm – 380 mm)With Diff. of 20 mm	EXCEL DNAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	289.25	G.T.N./ EXCEL SET	As per Annexure C	-	Excel nail YIPL/NL/ 03

	9 mm x (280 mm – 380 mm) With Diff. of 20	EXCEL DNAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	289.25	G.T.N./ EXCEL SET	As per Annexure C	-	Excel nail YIPL/NL/03
	10 mm x (280 mm – 380 mm) With Diff. of 20 mm	EXCEL DNAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	289.25	G.T.N./ EXCEL SET	As per Annexure C	-	Excel nail YIPL/NL/03
CISTA NAILS - LOCKING SYSTEM	8 mm x (280 mm – 380 mm) With Diff. of 20 mm	CISTA NAIL - S.S.	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs) & 3.9 mm (3 Pcs)	-	800.01	Spare	Cista Jig, Cista Bolt, Cista Attachment bolt & Cista Sleeve	-	Cista Nail YIPL/NL/02
	9 mm x (280 mm – 380 mm) With Diff. of 20 mm	CISTA NAIL - S.S	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (8 Pcs)	-	800.01	Spare	Cista Jig, Cista Bolt, Cista Attachment bolt & Cista Sleeve	-	Cista nail YIPL/NL/02
	10 mm x (280 mm – 380 mm) With Diff. of 20 mm	CISTA NAIL - S.S	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (8 Pcs)	-	800.01	Spare	Cista Jig, Cista Bolt, Cista Attachment bolt & Cista Sleeve	-	Cista Nail YIPL/NL/02

FEMUR NAIL LOCKING SYSTEM	9 mm x (360 mm – 440 mm)& Diff. of 10 mm	INTERLOCKING NAIL FOR FEMUR	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.31	Tibia/Femur Set	As per Annexure A	-	Interlocking nail for femur YIPL/NL /05
	10 mm x (360 mm – 440 mm)With Diff. of 10 mm	INTERLOCKING NAIL FOR FEMUR	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.31	Tibia/Femur Set	As per Annexure A	-	Interlocking nail for femur YIPL/NL /05
	11 mm x (360 mm – 440 mm)With Diff. of 10 mm	INTERLOCKING NAIL FOR FEMUR	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.31	Tibia/Femur Set	As per Annexure A	-	Interlocking nail for femur YIPL/NL /05
	12 mm x (360 mm – 440 mm)With Diff. of 10 mm	INTERLOCKING NAIL FOR FEMUR	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.31	Tibia/Femur Set	As per Annexure A	-	Interlocking nail for femur YIPL/NL /05
R- FEMUR NAIL LOCKING SYSTEM	9 mm x (360 mm – 420 mm)With Diff. of 20 mm	R- FEMUR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.31	Tibia/Femur Set	As per Annexure A	-	R femur nail YIPL/ST/ 03

	10 mm x (360 mm – 420 mm)With Diff. of 20 mm	R- FEMUR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.31	Tibia/Femur Set	As per Annexure A	-	R femur nail YIPL/ST/03
	11 mm x (360 mm – 420 mm)With Diff. of 20 mm	R- FEMUR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.31	Tibia/Femur Set	As per Annexure A	-	R femur nail YIPL/ST/03
	12 mm x (360 mm – 420 mm)With Diff. of 20 mm	R- FEMUR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.31	Tibia/Femur Set	As per Annexure A	-	R femur nail YIPL/ST/03
SUPRACONDYLAR NAIL LOCKING SYSTEM	9 mm x (150 mm – 350 mm) With Diff. of 50 mm	SUPRACONDYLAR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (3 Pcs) & 6.5 CANC. SCREW (2 Pcs)	-	291.12	SUPRA SET	As per Annexure D	-	Supracondylar nail YIPL/NL/08
	10 mm x (150 mm – 350 mm)With Diff. of 50 mm	SUPRACONDYLAR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (3 Pcs) & 6.5 CANC. SCREW (2 Pcs)	-	291.12	SUPRA SET	As per Annexure D	-	Supracondylar nail YIPL/NL/08

	11 mm x (150 mm – 350 mm) With Diff. of 50 mm	SUPRACON DYLAR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt 4.9 mm (3 Pcs) & 6.5 CANC. SCREW (2 Pcs)	-	291.12	SUPRA SET	As per Annexure D	-	Supracon dylar nail YIPL/NL /08
	12 mm x (150 mm – 350 mm) With Diff. of 50 mm	SUPRACON DYLAR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt 4.9 mm (3 Pcs) & 6.5 CANC. SCREW (2 Pcs)	-	291.12	SUPRA SET	As per Annexure D	-	Supracon dylar nail YIPL/NL /08
PFN - SHORT-130° LOCKING SYSTEM	9 mm x 250 mm	PROXIMAL FEMUR NAIL - SHORT-130*	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2
	10 mm x 250 mm	PROXIMAL FEMUR NAIL - SHORT-130*	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2

	11 mm x 250 mm	PROXIMAL FEMUR NAIL - SHORT-130*	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2
	12 mm x 250 mm	PROXIMAL FEMUR NAIL - SHORT-130*	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2
PFN - SHORT-135 ⁰ LOCKING SYSTEM	9 mm x 250 mm	PROXIMAL FEMUR NAIL - SHORT-135*	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2
	10 mm x 250 mm	PROXIMAL FEMUR NAIL - SHORT-135*	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2

	11 mm x 250 mm	PROXIMAL FEMUR NAIL - SHORT-135*	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2
	12 mm x 250 mm	PROXIMAL FEMUR NAIL - SHORT-135*	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2
TRAFFON NAIL - 130° LOCKING SYSTEM	9 mm x 180 mm	TRAFFON NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	Proximal femur nail YIPL/ST/0 2
	10 mm x 180 mm	TRAFFON NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	Proximal femur nail YIPL/ST/0 2

	11 mm x 180 mm	TRAFFON NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	Proximal femur nail YIPL/ST/O 2
	12 mm x 180 mm	TRAFFON NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	Proximal femur nail YIPL/ST/O 2
TRAFFON NAIL - 135° LOCKING SYSTEM	9 mm x 180 mm	TRAFFON NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	YIPL/ST/O 2
	10 mm x 180 mm	TRAFFON NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	Proximal femur nail YIPL/ST/O 2

	11 mm x 180 mm	TRAFFON NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	Proximal femur nail YIPL/ST/0 2
	12 mm x 180 mm	TRAFFON NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	Proximal femur nail YIPL/ST/0 2
LONG P.F.NAIL - 130° - LOCKING SYSTEM	9 mm x (340 mm - 440 mm) With Diff. of 20 mm	LONG P.F.NAIL - 130° - RIGHT	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2
	10 mm x (340 mm - 440 mm) With Diff. of 20 mm	LONG P.F.NAIL - 130° - RIGHT	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/0 2

	11 mm x (340 mm – 440 mm) With Diff. of 20 mm	LONG P.F.NAIL - 130° - RIGHT	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/O 2
LONG P.F.NAIL - 135* - LOCKING SYSTEM	9 mm x (340 mm – 440 mm) With Diff. of 20 mm	LONG P.F.NAIL - 130° - RIGH	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/O 2
	10 mm x (340 mm – 440 mm) With Diff. of 20 mm	LONG P.F.NAIL - 135* - RIGHT	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/O 2
	11 mm x (340 mm – 440 mm) With Diff. of 20 mm	LONG P.F.NAIL - 135* - RIGHT	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt 4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/O 2

	12 mm x (340 mm – 440 mm) With Diff. of 20 mm	LONG P.F.NAIL - 130° - LEFT	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/O 2
LONG P.F.NAIL - 135* - LOCKING SYSTEM	9 mm x (340 mm – 440 mm) With Diff. of 20 mm	LONG P.F.NAIL - 135* - LEFT	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/O 2
	10 mm x (340 mm – 440 mm) With Diff. of 20	LONG P.F.NAIL - 135* - LEFT	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/O 2
	11 mm x (340 mm – 440 mm) With Diff. of 20 mm	LONG P.F.NAIL - 135* - LEFT	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur nail YIPL/ST/O 2

HUMERUS -1 D.H. NAIL LOCKING SYSTEM	6 mm x (200 mm – 300 mm) With Diff. of 10 mm	INTERLOCKING NAIL FOR HUMERUS - 1 D.H.	-	-	-	Locking Bolt 2.9 mm (2 Pcs)	-	288.11	HUMERUS SET	As per Annexure F	-	Humerus Nail YIPL/NL /04
	7 mm x (200 mm – 300 mm) With Diff. of 10 mm	INTERLOCKING NAIL FOR HUMERUS - 1 D.H.	-	-	Guide Wire 2.0 x 600 (2 pcs)	Locking Bolt 3.9 mm (2 Pcs)	-	288.11	HUMERUS SET	As per Annexure F	-	Humerus Nail YIPL/NL /04
	8 mm x (200 mm – 300 mm) v Diff. of 10 mm	INTERLOCKING NAIL FOR HUMERUS - 1 D.H.	-	-	Guide Wire 2.0 x 600 (2 pcs)	Locking Bolt 3.9 mm (2 Pcs)	-	288.11	HUMERUS SET	As per Annexure F	-	Humerus Nail YIPL/NL /04
PROXIMAL HUMERUS NAIL - LOCKING SYSTEM	7 mm x (150 mm – 300 mm) With Diff. of 50 mm	PROXIMAL HUMERUS NAIL - 2DH CANN	-	-	Guide Wire 2.0 x 600 (2 pcs)	Locking Bolt 3.9 mm (3 Pcs)	-	288.11	HUMERUS SET	As per Annexure F	-	Humerus Nail YIPL/NL /04
	8 mm x (150 mm – 300 mm) With Diff. of 50 mm	PROXIMAL HUMERUS NAIL - 2DH CANN	-	-	Guide Wire 2.0 x 600 (2 pcs)	Locking Bolt 3.9 mm (3 Pcs)	-	288.11	HUMERUS SET	As per Annexure F	-	Humerus Nail YIPL/NL /04

	9 mm x (150 mm – 300 mm) With Diff. of 50 mm	PROXIMAL HUMERUS NAIL - 2DH CANN	-	-	Guide Wire 2.0 x 600 (2 pcs)	Locking Bolt4.9 mm (3 Pcs)	-	288.11	HUMERUS SET	As per Annexure F	-	Humerus Nail YIPL/NL /04
R- HUMERUS NAIL - LOCKING SYSTEM	6 mm x 205, 295, 300 mm, From (220 mm – 280 mm) With Diff. of 10 mm	R HUMERUS NAIL – SOLID	-	-	-	Locking Bolt3.4 mm (6 Pcs)	-	288.60	R- HUMERUS SET	As per Annexure G	-	Humerus Nail YIPL/NL /04
	6.7 mm x 205, 295, 300 mm, From (220 mm – 280 mm) With Diff. of 10 mm	R HUMERUS NAIL – SOLID	-	-	-	Locking Bolt3.4 mm (6 Pcs)	-	288.60	R- HUMERUS SET	As per Annexure G	-	Humerus Nail YIPL/NL /04
	7.5 mm x 205, 295, 300 mm, From (220 mm – 280 mm) With Diff. of 10 mm	R HUMERUS NAIL – SOLID	-	-	-	Locking Bolt3.4 mm (6 Pcs)	-	288.60	R- HUMERUS SET	As per Annexure G	-	Humerus Nail YIPL/NL /04

ELASTIC NAIL LOCKING SYSTEM	2 mm X 350 mm	ELASTIC NAIL	-	-	-	-	-	500.10	ELASTIC SET	As per Annexure H	-	Elastic Nail YIPL/NL /01
	2.5 mm X 400 mm	ELASTIC NAIL	-	-	-	-	-	500.10	ELASTIC SET	As per Annexure H	-	Elastic Nail YIPL/NL /01
	3.0 mm X 400 mm	ELASTIC NAIL	-	-	-	-	-	500.10	ELASTIC SET	As per Annexure H	-	Elastic Nail YIPL/NL /01
	3.5 mm X 400 mm	ELASTIC NAIL	-	-	-	-	-	500.10	ELASTIC SET	As per Annexure H	-	Elastic Nail YIPL/NL /01
	4.0 mm X 450 mm	ELASTIC NAIL	-	-	-	-	-	500.10	ELASTIC SET	As per Annexure H	-	Elastic Nail YIPL/NL /01
DISTAL FEMUR NAIL LOCKING SYSTEM	10 mm x (150 mm – 250 mm) With Diff. of 50 mm	DISTAL FEMUR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt 4.9 mm (4 Pcs) & 6.5 CANC. SCREW (2 Pcs)	-	291.12	SUPRA SET	As per Annexure D	-	Supracondylar nail YIPL/NL /08
	10 mm x (300 mm – 420 mm) With Diff. of 20 mm	DISTAL FEMUR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt 4.9 mm (4 Pcs) & 6.5 CANC. SCREW (2 Pcs)	-	291.12	SUPRA SET	As per Annexure D	-	Supracondylar nail YIPL/NL /08

	11 mm x (150 mm – 250 mm) With Diff. of 50 mm	DISTAL FEMUR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs) & 6.5 CANC. SCREW (2 Pcs)	-	291.12	SUPRA SET	As per Annexure D	-	Supracon dylar nail YIPL/NL /08
	11 mm x (300 mm – 420 mm) With Diff. of 20 mm	DISTAL FEMUR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs) & 6.5 CANC. SCREW (2 Pcs)	-	291.12	SUPRA SET	As per Annexure D	-	Supracon dylar nail YIPL/NL /08
	12 mm x (150 mm – 350 mm) With Diff. of 50 mm	DISTAL FEMUR NAIL	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs) & 6.5 CANC. SCREW (2 Pcs)	-	291.12	SUPRA SET	As per Annexure D	-	Supracon dylar nail YIPL/NL /08
TRAFFON II (SMALL) LOCKING SYSTEM	9 mm X 180 mm	S.S. TRAFFON II (SMALL) – 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
	10 mm X 180 mm	S.S. TRAFFON II (SMALL) – 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
	11 mm X 180 mm	S.S. TRAFFON II (SMALL) – 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09

	12 mm X 180 mm	S.S. TRAFFON II (SMALL) – 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
S.S. TRAFFON II (REGULA R) LOCKING SYSTEM	9 mm X 240 mm	S.S. TRAFFON II (REGULAR) – 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
	10 mm X 240 mm	S.S. TRAFFON II (REGULAR) – 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
	11 mm X 240 mm	S.S. TRAFFON II (REGULAR) – 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
	12 mm X 240 mm	S.S. TRAFFON II (REGULAR) – 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
LONG TRAFFON II NAIL- RIGHT LOCKING SYSTEM	9 mm x (300 mm – 420 mm) With Diff. of 20 mm	S.S. LONG TRAFFON II NAIL - 130° - RIGHT	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09

	10 mm x (300 mm – 420 mm) With Diff. of 20 mm	S.S. LONG TRAFFON II NAIL - 130° - RIGHT	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFON II YIPL/ST/08
	11 mm x (300 mm – 420 mm) With Diff. of 20 mm	S.S. LONG TRAFFON II NAIL - 130° - RIGHT	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFON II YIPL/NL/09
LONG TRAFFON II NAIL - LEFT LOCKING SYSTEM	9 mm x (300 mm – 420 mm) With Diff. of 20 mm	S.S. LONG TRAFFON II NAIL - 130° - LEFT	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFON II YIPL/NL/09
	10 mm x (300 mm – 420 mm) With Diff. of 20 mm	S.S. LONG TRAFFON II NAIL - 130° - LEFT	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFON II YIPL/NL/09
	11 mm x (300 mm – 420 mm) With Diff. of 20 mm	S.S. LONG TRAFFON II NAIL - 130° - LEFT	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFON II YIPL/NL/09

SQUARE NAIL LOCKING SYSTEM	2.0 mm X 300 mm	SQUARE NAIL FOR ULNA & RADIUS - UPTO - 300MM	-	-	-	-	-	-	-	-	-	Square Nail YIPL/NL /07
	2.5 mm X 300 mm	SQUARE NAIL FOR ULNA & RADIUS - UPTO - 300MM	-	-	-	-	-	-	-	-	-	Square Nail YIPL/NL /07
	3.0 mm X 300 mm	SQUARE NAIL FOR ULNA & RADIUS - UPTO - 300MM	-	-	-	-	-	-	-	-	-	Square Nail YIPL/NL /07
	3.5 mm X 300 mm	SQUARE NAIL FOR ULNA & RADIUS - UPTO - 300MM	-	-	-	-	-	-	-	-	-	Square Nail YIPL/NL /07
	4.0 mm X 300 mm	SQUARE NAIL FOR ULNA & RADIUS - UPTO - 300MM	-	-	-	-	-	-	-	-	-	Square Nail YIPL/NL /07

K ' NAILFEMUR LOCKING SYSTEM	9 mm X 440 mm	K ' NAIL FOR FEMUR	-	-	-	-	-	-	-	-	-	-
	10 mm X 440 mm	K ' NAIL FOR FEMUR	-	-	-	-	-	-	-	-	-	-
	11 mm X 440 mm	K ' NAIL FOR FEMUR	-	-	-	-	-	-	-	-	-	-
	12 mm X 440 mm	K ' NAIL FOR FEMUR	-	-	-	-	-	-	-	-	-	-
V ' NAIL TIBIA LOCKING SYSTEM	8 mm X 380 mm	V ' NAIL FOR TIBIA	-	-	-	-	-	-	-	-	-	-
	9 mm X 380 mm	V ' NAIL FOR TIBIA	-	-	-	-	-	-	-	-	-	-
	10 mm X 380 mm	V ' NAIL FOR TIBIA	-	-	-	-	-	-	-	-	-	-
CISTA NAILS - TITANIUM	8 mm x (280 mm – 380 mm) With	CISTA NAILS - TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2	Locking Bolt4.9 mm (5 Pcs) & 3.9 mm	-	800.01	Spare	Cista Jig, Cista Bolt, Cista Attachment	-	Cista nail YIPL/NL /02

LOCKING SYSTEM	Diff. of 20 mm				pcs)	(3 Pcs)				bolt & Cista Sleeve		
	9 mm x (280 mm – 380 mm) With Diff. of 20 mm	CISTA NAILS - TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt 4.9 mm (5 Pcs) & 3.9 mm (3 Pcs)	-	800.01	Spare	Cista Jig, Cista Bolt, Cista Attachment bolt & Cista Sleeve	-	Cista nail YIPL/NL /02
	10 mm x (280 mm – 380 mm) With Diff. of 20 mm	CISTA NAILS - TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt 4.9 mm (5 Pcs) & 3.9 mm (3 Pcs)	-	800.01	Spare	Cista Jig, Cista Bolt, Cista Attachment bolt & Cista Sleeve	-	Cista nail YIPL/NL /02
R- TIBIA NAIL TITANIUM LOCKING SYSTEM	8 mm x (280 mm – 380 mm) With Diff. of 20 mm	R- TIBIA NAIL TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt 3.9 mm (5 Pcs)	-	287.49	R-TIBIA SET	As per Annexure J	-	Expert nail YIPL/NL /03
	9 mm x (280 mm – 380 mm) With Diff. of 20 mm	R- TIBIA NAIL TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt 4.9 mm (5 Pcs)	-	287.49	R-TIBIA SET	As per Annexure J	-	Expert nail YIPL/NL /03
	10 mm x (280 mm – 380 mm) With Diff. of 20 mm	R- TIBIA NAIL TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt 4.9 mm (5 Pcs)	-	287.49	R-TIBIA SET	As per Annexure J	-	Expert nail YIPL/NL /03

	20 mm											
R- FEMUR NAIL TITANIUM LOCKING SYSTEM	9 mm x (360 mm – 440 mm) With Diff. of 20 mm	R- FEMUR NAIL TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.48	R-FEMUR SET	As per Annexure A	-	Expert nail YIPL/NL /03
	10 mm x (360 mm – 440 mm) With Diff. of 20 mm	R- FEMUR NAIL TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.48	R-FEMUR SET	As per Annexure A	-	Expert nail YIPL/NL /03
	11 mm x (360 mm – 440 mm) With Diff. of 20 mm	R- FEMUR NAIL TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.48	R-FEMUR SET	As per Annexure A	-	Expert nail YIPL/NL /03
	12 mm x (360 mm – 440 mm) With Diff. of 20 mm	R- FEMUR NAIL TITANIUM	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (4 Pcs)	-	287.48	R-FEMUR SET	As per Annexure A	-	Expert nail YIPL/NL /03
TITANIUM FACT GT - RIGHT	9 mm x (340 mm – 440 mm) With	TITANIUM FACT GT - RIGHT	-	-	Guide Wire 2.5 x 900 mm(2	Locking Bolt4.9 mm (5 Pcs), 6.4 Proximal	-	700.25	FACT GT SET	As per Annexure K	-	Proximal femoral Nail YIPL/ST/

LOCKING SYSTEM	Diff. of 20 mm				pcs)	Bolt (2 Pcs)						02
	10 mm x (340 mm – 440 mm) With Diff. of 20 mm	TITANIUM FACT GT - RIGHT	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs), 6.4 Proximal Bolt (2 Pcs)	-	700.25	FACT GT SET	As per Annexure K	-	Proximal femoral Nail YIPL/ST/02
	11 mm x (340 mm – 440 mm) With Diff. of 20 mm	TITANIUM FACT GT - RIGHT	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs), 6.4 Proximal Bolt (2 Pcs)	-	700.25	FACT GT SET	As per Annexure K	-	Proximal femoral Nail YIPL/ST/02
	12 mm x (340 mm – 440 mm) With Diff. of 20 mm	TITANIUM FACT GT - RIGHT	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs), 6.4 Proximal Bolt (2 Pcs)	-	700.25	FACT GT SET	As per Annexure K	-	Proximal femoral Nail YIPL/ST/02
TITANIUM FACT GT - LEFT LOCKING SYSTEM	9 mm x (340 mm – 440 mm) With Diff. of 20 mm	TITANIUM FACT GT - LEFT	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs), 6.4 Proximal Bolt (2 Pcs)	-	700.25	FACT GT SET	As per Annexure K	-	Proximal femoral Nail YIPL/ST/02
	10 mm x (340 mm – 440 mm) With Diff. of	TITANIUM FACT GT - LEFT	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs), 6.4 Proximal Bolt	-	700.25	FACT GT SET	As per Annexure K	-	Proximal femoral Nail YIPL/ST/02

	20 mm					(2 Pcs)						
	11 mm x (340 mm – 440 mm) With Diff. of 20 mm	TITANIUM FACT GT - LEFT	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs), 6.4 Proximal Bolt (2 Pcs)	-	700.25	FACT GT SET	As per Annexure K	-	Proximal femoral Nail YIPL/ST/02
	12 mm x (340 mm – 440 mm) With Diff. of 20 mm	TITANIUM FACT GT - LEFT	-	-	Guide Wire 2.5 x 900 mm(2 pcs)	Locking Bolt4.9 mm (5 Pcs), 6.4 Proximal Bolt (2 Pcs)	-	700.25	FACT GT SET	As per Annexure K	-	Proximal femoral Nail YIPL/ST/02
TITANIUM PFN NAIL - LOCKING SYSTEM	9 mm X 250 mm	TITANIUM PFN NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur Nail YIPL/ST/02
	10 mm X 250 mm	TITANIUM PFN NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur Nail YIPL/ST/02

	11 mm X 250 mm	TITANIUM PFN NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur Nail YIPL/ST/ 02
	12 mm X 250 mm	TITANIUM PFN NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur Nail YIPL/ST/ 02
TITANIUM PFN NAIL - LOCKING SYSTEM	9 mm X 250 mm	TITANIUM PFN NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur Nail YIPL/ST/ 02
	10 mm X 250 mm	TITANIUM PFN NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur Nail YIPL/ST/ 02

	11 mm X 250 mm	TITANIUM PFN NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur Nail YIPL/ST/ 02
	12 mm X 250 mm	TITANIUM PFN NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.35	PFN SET	As per Annexure E	-	Proximal femur Nail YIPL/ST/ 02
TITANIUM TRAFFON NAIL - LOCKING SYSTEM	9 mm X 180 mm	TITANIUM TRAFFON NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	TRAFFR ON II YIPL/NL /09
	10 mm X 180 mm	TITANIUM TRAFFON NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	TRAFFR ON II YIPL/NL /09

	11 mm X 180 mm	TITANIUM TRAFFON NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	TRAFFR ON II YIPL/NL /09
	12 mm X 180 mm	TITANIUM TRAFFON NAIL - 135°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	TRAFFR ON II YIPL/NL /09
TITANIUM TRAFFON NAIL - LOCKING SYSTEM	9 mm X 180 mm	TITANIUM TRAFFON NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	TRAFFR ON II YIPL/NL /09
	10 mm X 180 mm	TITANIUM TRAFFON NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	TRAFFR ON II YIPL/NL /09

	11 mm X 180 mm	TITANIUM TRAFFON NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	TRAFFR ON II YIPL/NL /09
	12 mm X 180 mm	TITANIUM TRAFFON NAIL - 130°	-	-	Guide Wire 2.5 x 400 (2 pcs) & 2.0 x 400 mm(2 pcs)	Locking Bolt4.9 mm (2 Pcs), 6.4 Proximal Bolt (1 Pcs) & 8 mm PFN Bolt (1 pcs)	-	293.36	TFN SET	As per Annexure F	-	TRAFFR ON II YIPL/NL /09
POLARIS NAIL LOCKING SYSTEM	6 mm X 150 mm	POLARIS NAIL	-	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt2.9 mm (6 Pcs),	-	288.50	POLARIS SET	As per Annexure L	-	Humerus nail YIPL/NL /04
	7 mm X 150 mm	POLARIS NAIL	-	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt3.9 mm (6 Pcs),	-	288.50	POLARIS SET	As per Annexure L	-	Humerus nail YIPL/NL /04
	8 mm X 150 mm	POLARIS NAIL	-	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt3.9 mm (6 Pcs),	-	288.50	POLARIS SET	As per Annexure L	-	Humerus nail YIPL/NL /04

	9 mm X 150 mm	POLARIS NAIL	-	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (6 Pcs),	-	288.50	POLARIS SET	As per Annexure L	-	Humerus nail YIPL/NL /04
POLARIS PLUS NAIL LOCKING SYSTEM	6 mm x (200 mm – 300 mm) With Diff. of 20 mm	POLARIS PLUS NAIL	-	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt2.9 mm (6 Pcs),	-	288.50	POLARIS SET	As per Annexure L	-	Humerus nail YIPL/NL /04
	7 mm x (200 mm – 300 mm) With Diff. of 20 mm	POLARIS PLUS NAIL	-	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt3.9 mm (6 Pcs),	-	288.50	POLARIS SET	As per Annexure L	-	Humerus nail YIPL/NL /04
	8 mm x (200 mm – 300 mm) With Diff. of 20 mm	POLARIS PLUS NAIL	-	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt3.9 mm (6 Pcs),	-	288.50	POLARIS SET	As per Annexure L	-	Humerus nail YIPL/NL /04
	9 mm x (200 mm – 300 mm) With Diff. of 20 mm	POLARIS PLUS NAIL	-	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (6 Pcs),	-	288.50	POLARIS SET	As per Annexure L	-	Humerus nail YIPL/NL /04
TITANIUM TRAFFON II	9 mm x 180 mm	TITANIUM TRAFFON II (SMALL) – 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFON II YIPL/NL /09

(SMALL) LOCKING SYSTEM					pcs)							
	10 mm x 180 mm	TITANIUM TRAFFON II (SMALL) - 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
	11 mm x 180 mm	TITANIUM TRAFFON II (SMALL) - 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
	12 mm x 180 mm	TITANIUM TRAFFON II (SMALL) - 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
TITANIUM TRAFFON II (REGULA R) LOCKING SYSTEM	9 mm x 240 mm	TITANIUM TRAFFON II (REGULAR) - 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09
	10 mm x 240 mm	TITANIUM TRAFFON II (REGULAR) - 130 ⁰	Proxi mal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFR ON II YIPL/NL /09

	11 mm x 240 mm	TITANIUM TRAFFON II (REGULAR) - 130°	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFRON II YIPL/NL /09
	12 mm x 240 mm	TITANIUM TRAFFON II (REGULAR) - 130°	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFRON II YIPL/NL /09
TITANIUM LONG TRAFFON II NAIL - RIGHT LOCKING SYSTEM	9 mm x 300 mm, (340 mm – 420 mm) With Diff. of 20 mm	TITANIUM LONG TRAFFON II NAIL - 130° - RIGHT	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFRON II YIPL/NL /09
	10 mm x 300 mm, (340 mm – 420 mm) With Diff. of 20 mm	TITANIUM LONG TRAFFON II NAIL - 130° - RIGHT	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFRON II YIPL/NL /09
	11 mm x 300 mm, (340 mm – 420 mm) With Diff. of 20 mm	TITANIUM LONG TRAFFON II NAIL - 130° - RIGHT	Proximal Blade (1 pcs)	-	Guide Wire 2.5 x 600 mm(2 pcs)	Locking Bolt4.9 mm (1 Pcs)	-	501.65	TRAFFON II SET	As per Annexure I	-	TRAFFRON II YIPL/NL /09

